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- **Cleaning**
- **Disinfecting**
- **Storing**
- **Use at baths and steam baths**
- **Use at swimming pools**

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## Care instructions

Dear Customer

Congratulations on your purchase of a high-quality AIREX brand product! To enjoy the benefits of our CE-conform AIREX speciality foam products for as long as possible, please note these care instructions.

### Cleaning

Dirty mats can be cleaned with a soft brush and hand-hot soapy water or with a solution of household detergent (washing-up liquid). Although washing in a washing machine is possible, it is not recommended for practical reasons. If unavoidable, select a gentle wash program and do not exceed a temperature of 40°C. Standard laundry detergent should be used. Benzine, alcohol and other standard stain removers should never be used for cleaning under any circumstances.

### Disinfecting

Generally speaking, AIREX exercise mats do not require disinfecting thanks to their antibacterial finish. However, if disinfection is required for special hygienic reasons, an aqueous solution of any of the standard disinfectants and detergents found in households, hospitals and institutes should be used (Merfen, Descogen, Kohrsolin, Kohrsolin FF, Terralin). Disinfection with hot steam is not permissible.

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### Storing

Individual mats can be simply rolled up and stored with the straps which can be supplied for this purpose. For the simultaneous storage of several mats, draping over a mat stand or flat storage are both ideal. If left unused for longer periods of time, stacked mats should preferably be stored in dry, unheated rooms.

### Use at baths and steam baths

In order to ensure AIREX mats enjoy a long service life, they should be rinsed thoroughly at regular intervals and air dried when used in warm and extremely damp surroundings such as public and private baths, spas or similar.

### Use at swimming pools

AIREX mats are resistant to standard water treatment chemicals and bactericides (chlorine, hypochlorite, ozone, etc.). Rinsing the mats after use with fresh water and air drying are recommended (also for hygienic reasons). In principle, the mats are resistant to vegetable oils and greases. However, some sun protection oils and emulsions may slightly affect the surfaces and thus reduce their service life. The use of a textile cover, e.g. a bath towel, is recommended in such cases.

### AIREX Gymnastic and Rehabilitation mats are:

- safe
- comfortable and body-friendly
- hygienic
- user-friendly
- ecologically manufactured

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